

QEL, Spices Board


*In-house Training – 1:
Analysis of Pesticide Residues*

Lecture - 2

*Introduction to instrumentation:
GC-MS/MS*

10-Jul-23

Dr. Ramesh BN
Spices Board



Moving on

From last lecture...

What we have **already** covered:

- ↳ Basics of chromatography and HPLC
- ↳ General principles of mass spectrometry (MS) and hyphenated techniques
- ↳ General principles of LC-MS/MS
 - ↳ Quadrupole theory
 - ↳ How mass filtration happens in MS: scanning (RF-DC ramp) and selection (RF-DC fixed)

What we **will** cover today:

- ↳ Review of basics operation of **GC**
- ↳ General principles of **GC-MS** and **GC-MS/MS**
- ↳ How GC-MS differs from **LC-MS**
- ↳ Using GC-MS and GC-MS/MS for residue analysis


Gas Chromatography

A review...

- ⇒ Chromatographic technique for **separation** of mixtures of **volatile compounds**.
 - ↳ Mobile Phase is a **gas**, He or N₂.
 - ↳ Mixtures are separated based on the difference in interaction of the components with the stationary phase.
 - ↳ Identification of unknowns is based on matching TR of standards
- ⇒ Using **conventional** detectors require relying on **retention time matching** for identification causes problems:
 - ↳ False **positives**, false **negatives**
- ⇒ **Mass spectrometry** is used as the detector to avoid these problems

Gas Chromatography

parts




The components

- ↳ Helium gas is the carrier
- ↳ Moisture, oxygen and hydrocarbon traps
- ↳ The column is a capillary tube of fused silica, with the liquid stationary phase coated inside
- ↳ The electronic flow controller (EFC) manages the flow of carrier gas through the column
- ↳ The oven controls the temperature of the column
- ↳ The injector vaporizes the sample and mixes it with the carrier gas
- ↳ Detector detects the compounds separated in the column


MS as a detector for GC

Sensitivity and range compared to conventional detectors



The carrier gas

Things to remember




- ⇒ The carrier gas cylinder should be replaced when the inner pressure falls below 200 psi
 - ↳ Will affect MS performance if pressure is lower

GC-MS and –MS/MS

§ Components and operation

The interface

Hyphenation to GC



- Column outlet is placed directly opposite to the MS Inlet
- MS vacuum can handle around **1 ml/min carrier flow**
 - A flow-restrictor is used
 - High GC flow – excess is vented out
 - Low GC flow – He from the reservoir is pulled in
- Best for capillary columns

How GC-MS/MS works

General principles

⇒ Very similar to LC-MS/MS, which we have already covered

- ↳ MS can 'see' only ions, so **ionization** is the first step in MS
- ↳ In GC-MS/MS,
 - ↳ ionization happens in **vacuum** (unlike LC-MS)
 - ↳ ionization is '**hard**' (unlike LC-MS)
- ↳ Other operations are similar to LC-MS/MS, which has been covered already:
 - ↳ **Quadrupole** theory
 - ↳ Mass **scanning** / Mass **filtration** (RF/DC)
 - ↳ Collision induced dissociation
 - ↳ Detection by electron multiplier

The Vacuum system

In GC-MS


Like in LC-MS, in GC-MS also the vacuum is achieved in two stages:

- ↳ The **roughing pump** (outside the instrument): from **760 Torr** to **10^{-3} Torr**
- ↳ The **turbo molecular pump** (inside the instrument): from **10^{-3} Torr** to **10^{-5} Torr**



GC-MS/MS




Hyphenation to GC



In GC MS/MS, ionization happens in vacuum
Ionization is hard: Molecular ion mostly doesn't survive
Standard spectrum library is possible in GC MS mode

Ionization

In GC-MS



Mass spectrum
Which can be
library-matched

One or more of these ions
can become **parent ions**
for MS/MS analysis

GC-MS


Identification by library matching (1)

⇒ GC-MS ionization is **hard**

- ↳ The parent ion, M^+ , most often does survive
- ↳ It breaks down into **several smaller ions**
- ↳ These smaller ions can be **scanned** (one by one, from lowest to highest), and detected
- ↳ A **mass spectrum** is produced

⇒ Conditions are **standard** (unless we choose to change them):

- ↳ Ionization is at **70 eV**
- ↳ Vacuum is set at **10^{-5} torr**
- ↳ Under standard conditions **spectra produced are also standard**
- ↳ Can be matched with standard **libraries**
 - ↳ **NIST**, **Wily** Pesticide Library etc.




If **molecular ion** survives, it will show up as the ion with **highest mass**

NIST
National Institute of
Standards and Technology
U.S. Department of Commerce

WILEY REGISTRY®
12th Edition/NIST 2020
Mass Spectral Library

GC-MS

Identification by library matching (2)




GC-MS is a very useful technique for identification (unlike LC-MS)

GC-MS/MS


Quantification by MRM mode(1)

- ⇒ The **base peak** (or any other intense ion) is selected as the **parent ion** from the mass spectrum
- ⇒ Subjected to collision induced dissociation
 - ↳ In the **second quadrupole**, N_2 is used as **collision gas**
 - ↳ The parent ion fragments into **daughter ions**
 - ↳ A daughter ion is selected in the **third quadrupole**, followed by detection
- ⇒ Multiple reaction monitoring (**MRM**)




GC-MS/MS

Quantification by MRM mode



“Transition”



The MRM table


Name	Transition	Dwell time
Compound A	$P_A \rightarrow D_1$	100 ms
Compound A	$P_A \rightarrow D_2$	100 ms
Compound B	$P_B \rightarrow D_1$	100 ms
Compound B	$P_B \rightarrow D_2$	100 ms
....		

In GC-MS, best transitions can vary from matrix to matrix

GC-MS/MS Quantification

Quantifier and qualifier

- In MRM, we choose 2 transitions per analyte
- Higher intensity transition is called the **quantifier**, and is used for quantification
- Lower intensity transition is called the **qualifier**
- Quantifier : Qualifier ratio should match between **standard and sample**



MRM time segmenting


Why?

Issue:


- GC peaks are typically **sharp**, usually only a few seconds wide
- When number of transitions in the MRM table is high, total **cycle time** will be high
- Number of **data points per peak** will be low
- Will affect **peak shapes** – a minimum of 5 points per peak needed

Solution:

- Segment the total run-time into **smaller periods**
- Within each period, set an MRM cycle
- Also called '**dynamic**' MRM




Number of data points high



Number of data points low


MRM time segmenting (2)

Also applicable in LC-MS/MS



Key difference between GC- and LC-MS/MS

Matrix effect



The calibration curve

Usually,

- Matrix **suppression** in LC-MS/MS
- Matrix **enhancement** in GC-MS/MS

We can't use solvent calibration curve for analytical work

Matrix-matched calibration curve required

Problem: availability of **blank matrices**


Practical aspects in GC-MS

§ Some specific topics

Autotune in GC-MS

Startup checks

- ⇒ When MS shows errors or low sensitivity, or if there has been an unforeseen shutdown, we perform an autotune
- ↳ involves adjusting **several mass spectrometer parameters** through the **infusion** of a tune compound, commonly perfluorotributylamine (**PFTBA**)
 - ↳ Usually the first step in troubleshooting




PFTBA

Autotune in GC-MS (2)

Startup checks

⇒ Autotune checks if:

- ↳ air/Water ratio is within limits (leaks and gas purity)
- ↳ background noise is low
- ↳ masses are correctly assigned (61, 219 and 502)
- ↳ the mass peak widths (PW50) are 0.55 ± 0.1 .
- ↳ the electron multiplier (EM) voltage is around 1500V. Higher values indicate saturation




Column trimming

Retention time changes

⇒ After large number of injections, head of the capillary column might become contaminated:

- ↳ Head of the column might need trimming
- ↳ After many trimmings, RT shift might occur
- ↳ Might need to adjust the segmentation in the method
- ↳ RT-locking calculators are available




Liner changing

Contamination issues

⇒ After large number of injections, if there is sudden drop in linearity, or tailing of peaks, liner might be the cause:

- ↳ Liners can get contaminated
- ↳ With spices, liner change might be needed within around 100 injections
- ↳ Regeneration and repacking of liners possible



More to learn..

Next session 17th June 2023

- ☑ (1) Introduction to instrumentation: LC-MS/MS, 03-Jul-23
- ☑ (2) Introduction to instrumentation: GC-MS/MS, 10-Jul-23
- 🕒 (3) Pesticide residue analysis - Introduction, 17-Jul-23 ←
- 🕒 (4) Advanced pesticide residue analysis, 24-Jul-23
- 🕒 (5) Method validation: requirements and practice, 31-Jul-23
- 🕒 (6) Introduction to measurement uncertainty calculation, 7-Aug-23

Dr Ramesh BN

Scientist C, Spices Board

WhatsApp : +91 8547662471

Mobile : +91 8547662471

Email (O) : rameshbabu.n@nic.in

Email (P) : rameshbn1@gmail.com

Course website : <http://ramesh-sb.github.io>

Thank you!
§ Questions?